

SOLID PLUS & SOLID PLUS HEAT LED AED PLACEMENT ASSISTANT

Place the defibrillator (AED) on the **platform**.
Use the **elastic band** to secure the defibrillator.

Cardiac Science

Powerheart G5

CU Medical Systems

I-Pad CU-SP1

I-Pad NF1200

Defibtech

Lifeline

Lifeline Pro

Lifeline View

HeartSine

Samaritan Pad 350P

Samaritan Pad 360P

Samaritan Pad 450P

Mediana

A10

A15

Mindray

BeneHeart D1

BeneHeart D1 Pro

Nihon Kohden

Cardiolife AED-2150

Cardiolife AED-3100

Philips

Heartstart FR3

Heartstart FRx

Heartstart OnSite/HS1

Physio control

CR2

Lifepak 1000

Lifepak CR Plus

Lifepak Express

Primedic

Heartsave AED

Heartsave AS

Heartsave One

Heartsave Pad

Progetti

Rescue Sam

Rescue Sam Pro

Radian

HR-501

HR-503

HR-503A

Schiller

Fred Easy

Fred Easy Life

Fred Easyport

Fred PA-1

Zoll

AED 3

AED Plus



SOLID PLUS & SOLID PLUS HEAT LED AED PLACEMENT ASSISTANT

Use **'both bolts'** (2x DIN screw M6) to hang the AED.
If needed, the AED can also be placed on a **'special bracket'**.



Cardiac Science

Powerheart G3 Plus
Powerheart G3 Pro
Powerheart G5

Cardiaid

CT0207RF
CT0207RS

CU Medical Systems

I-Pad CU-SP1
I-Pad NF1200

Defibtech

Lifeline
Lifeline Pro
Lifeline View

HeartSine

Samaritan Pad 350P
Samaritan Pad 360P
Samaritan Pad 450P

Mediana

A10
A15

Mindray

BeneHeart D1
BeneHeart D1 Pro

Nihon Kohden

Cardiolife AED-2150
Cardiolife AED-3100

Philips

Heartstart FR3
Heartstart FRx
Heartstart OnSite/HS1

Physio control

CR2
Lifepak 1000
Lifepak CR Plus
Lifepak Express

Primedic

Heartsave AED
Heartsave AS
Heartsave One
Heartsave Pad

Progetti

Rescue Sam
Rescue Sam Pro

Radian

HR-501
HR-503
HR-503A

Schiller

Fred Easy
Fred Easy Life
Fred Easyport
Fred PA-1

Zoll

AED 3
AED Plus

